





**CONTECH**<sup>®</sup>  
**Instruments Ltd.**  
 AN ISO 9001-2000 COMPANY  
*Where innovation never stops...*

## CARAT BALANCES

0.0005 ct to 1500 ct



- Bright Display.
- Standard bi-directional RS-232 interface.
- Multiple weighing units & functions.
- Rechargeable battery (optional).
- Overload protection.

Model	Capacity	Readability	Pan Size
CT-20C	100 ct	0.005 ct	100 mm
CA-34C	150 ct	0.0005 ct	100 mm
CA-64C	300 ct	0.0005 ct	100 mm
CB-50C	250 ct	0.005 ct	100 mm
CB-220DC	250 / 1100 ct	0.005 / 0.025 ct	100 mm
CB-330DC	250 / 1500 ct	0.005 / 0.050 ct	100 mm