




CONTECH[®]
Instruments Ltd.

AN ISO 9001-2000 COMPANY

Where innovation never stops...



CUSTOM BUILT
HYBRID CIRCUIT

CA SERIES HIGH CAPACITY PRECISION BALANCES



- Rugged design for industrial applications.
- Bright LED display.
- Multiple weighing units.
- Standard bi-directional RS-232 interface.
- Polyfunction : % weighing, counting, fill mode, GSM etc.

Model	Capacity	Readability	PanSize
CB-6001	6000g	0.1g	180*220 mm
CA-12001	12000g	0.1g	