




# HIGH PRECISION LABORATORY BALANCES

## HIGH PRECISION BALANCES

0.001g to 500g



Customized, attractive  
LCD display



Models	CA-123	CA-223	CA-323	CA-413	CA-503
Capacity	120g	220g	320g	410g	500g
Readability	0.001g	0.001g	0.001g	0.001g	0.001g
Repeatability (+/-)	0.001 g	0.001 g	0.001 g	0.001 g	0.002g
Linearity (+/-)	0.001 g	0.001 g	0.002 g	0.002 g	0.002g
Pan Size(mm)	100	100	100	100	100

The balance is supplied with Breeze shield as standard.