


CONTECH[®]
Instruments Ltd.
 AN ISO 9001-2000 COMPANY
 Where innovation never stops...

**CT SERIES
 JEWELLERY BALANCE**
 0.01g TO 10000g



JEWELLERY BALANCE

- Bright LED Display.
- Standard bi-directional RS-232 interface.
- Multiple weighing units & functions.
- Rechargeable battery (optional).
- Overload protection.

Model	Capacity	Readability	Pan Size
CT-300	300g	0.01g	Ø 115 mm
CT-600	600g	0.01g	
CT-3001	3000g	0.1g	Ø 150 mm
CT-6001	6000g	0.1g	210 x 240 mm
CT-10002	10000g	0.2g	210 x 240 mm
CT-15005	15000g	0.5g	210 x 240 mm