

## ACCESSORIES

### Magnetic Stirrer MS-01

Stirring speed	: 0~2300 rpm at no load.
Base Diameter	: 120mm
Maximum stirring capacity	: 1000ml
Power supply	: 6V DC Adaptor
Dimension	: 120 dia x 50mm height
Meter Weight	: 500 gms.



### Electrode Holder EH-01

Dimensions (mm)	: 115(L) x 35(W) x 28(H)
Weight	: 22 g
Suitable for Electrode	: Dia 13mm~16mm




### Stirr Bead SB-01

Dimensions	: 27(L) x 10(W) x 6(H)
Weight	: 5g



### Printer CMP-01

Printing Method	: Thermal
Character Type	: Alpha Numeric
Character per Line	: 30 max.
Communication	: RS 232
Baud Rate	: 9600
Printing Speed	: 90mm per sec.
Paper Width	: 57mm
Paper Cutting	: Manual
Power Supply	: 8v,3A DC
Dimension (mm)	: 135(L) x 195(D) x 100(H)
Paper Width	: 57mm
Printing Speed	: 90mm per sec.
Dimensions (mm)	: 135(L) x 195(D) x 100(H)
Net Weight	: 800g approx.



### Power Supply SMPS-04

Input	: 100~300 V AC, 47-63 Hz
Output	: 9VDC, 500mA
Dimensions with pin (mm)	: 50(L) x 70(H) x 70(D)
Dimensions without pin (mm)	: 50(L) x 70(H) x 37(D)
Weight	: 100g



### Battery Backup BBEC-01

Input	: 100~300 V AC, 47-63 Hz
Output	: 9VDC, 500mA
Battery	: 6V, 1.2 Ah
Weight	: 600g
Deminsions	: 150(L) x 65 (w) x 40(H)



### Calibration Solutions

Solution	Corresponding Electrode	Description	Accuracy
CSPH 4.0	pH	pH 4.0	± 0.1 pH
CSPH 7.0		pH 7.0	± 0.1 pH
CSPH 9.2		pH 9.2	± 0.1 pH
CSCN 1413	Conductivity	1413 µs/cm	±0.5%

