

MINI PRINTER

Printer having serial port can be connected to all Contech scales with RS-232 interface for printing weight data. Different applications can generate various reports which can be printed on the serial printer. GLP software with inbuilt RTC allows for the composition of GLP procedures for calibration test and calibration report with date and time on a printer connected to the serial port of the balances.



Thermal



Dot Matrix



301, Punit Indl. Premises, Turbhe Naka, Navi Mumbai - 400 705.
Tel. : +91 22 6139 3000 (100 lines)
E-mail : sales@contechindia.in
Website: www.contechindia.com

FACTORY / NABL LAB



Plot No. EL-221 TTC Indl. Area, MIDC (Electronic Zone),
Mhape, Navi Mumbai - 400 710.